

SKYLINER II – DRAWS YOU INTO A NEVER-BEFORE-SEEN WELLNESS SENSATION.













Enjoy your personal massage program according to your individual needs and expectations. The next time you use the chair, press the memory button and activate the automatic memo mode to enjoy your personal massage.

MULTIPLE POSITIONS MASSAGE

Innovative acupressure point's sensor allows you to enjoy body massage in all parts of your body. By combining a variety of innovative methods, you can relax and reduce stress.

5 LEVELS OF MASSAGE WIDTH ADJUSTABLE

According to different body type, the back massage curve and shoulder position are automatically detected, and the massage range is adjusted according to different parts of the human body. When massaged under tapping and finger pressure, the massage width can be adjusted freely in the range of 6-21 CM with 5 levels.

5 LEVELS OF MASSAGE SPEED ADJUSTABLE

According to different massage needs, the massage speed can be selected in 5 levels. You can freely adjust the massage speed, design your massage according to your wishes and combine different massage experiences with exquisite quality to feel completely at ease.

3-LEVEL AIR PRESSURE INTENSITY ADJUSTABLE

With a total of 38 airbags and air pressure intensity adjustable in three levels, full relaxation from the shoulders to the feet

LED TOUCHSCREEN

7 inches super touch screen, equipped with new UL operate system, making the massage progress clear to be seen. Users can control the massage program by one button



TECHNICAL DATA

Dimensions upright	80 x 148 x 120 cm
Dimensions reclining	80 x 180 x 105 cm
Weight	141 kg
Voltage	220-240V ~ 50/60Hz
Nominal power	220 W
Rated runtime	20 minutes
braintronics®	21 minutes
Certificates	(€@¥



Width of seating area 52 cm | Width of shoulder area (inside) 52 cm

AVAILABLE COLOURS





timetorela

V

SKYLINER MASSAGE CHAIR





Casada International GmbH Obermeiers Feld 3 33104 Paderborn | Germany www.casada.com

braintronics®

SPECIALS



BRAINTRONICS®

Braintronics® technology synchronises your brainwaves through audio stimulation. Combined with a specially developed massage programme your stress-ridden body will be able to relax. No matter whether you'd like to relax, dream or learn – anything's possible.

SHIATSI

SHIATSU Shiatsu (finger press) is a form of body therapy originating in Japan. The massage technique consists of soft, rhythmic, far-reaching stretches and rotations. The aim of Shiatsu is to simulate individual parts of the body and to mobilise the muscles.


3D MASSAGE 3D

In this new technique, the massage robot not only moves upwards and downwards, but also forwards and backwards.

KNEADING Mag

Skin and muscles are either clasped and kneaded between the thumb and the index finger or with both hands. This massage technique is used in particular to ease tension.

TAPPING





The techniques of kneading and tapping are combined in order to relieve tension and stimulate blood circulation at the same time.

ROLLING


Pleasant alternation between relaxation and releasing the muscles – also perfectly suited for a final massage, for relaxation and recoverv of the muscles.



Deeply penetrating Carbon infrared heat, blood vessels expand, the blood flow increases and the muscle tone is reduced. In

addition, it carries out a highly beneficial effect on nerve cells, which are responsible for reporting pain to the brain.

SOLE ROLLER



With this automatic program the massage heads are working very close to your body, they massage the whole back down to your buttocks. The S-shaped mechanic is constructed after the natural anatomy of the spine.



AIR COMPRESSION Rising and receding airbags have a pumping

motion which produces natural muscle reflex es. This is particularly beneficial in the calf and foot area.

ZERO SPACE

For the lying-position the massage chair moves forward and therefore requires very little distance to the wall behind. Space saving and comfortable!

7FRO GRAVITY Ø

You put yourself in the position of an astro-ZERO-G naut - that "weightless feeling"



AUTOMATIC Specifically developed automatic programmes for the whole body, which offer a big range of different massage techniques. ionally you can also turn on the heating function, if you'd like.

sage area. ADJUSTABLE FOOTREST The footrest can be extended, this way it is

BACK AREA

adaptable to different body sizes.


An individual applicable to any size mas-







AIR AREA Whole body airbag massage. Different are-Whole body airbag massage as are individually selectable.



Massage speed is adjustable.

The distance between the both massage heads can be set (narrow, medium and wide)



BIUFTOOTH Wireless data transmission.

HUMAN HANDS FEELING

Specially designed massage heads imitate the hand motion of a professional masseur. The maximum stretch space is 6,5 cm, and the maximum stretch angle is 41 degrees, so the mechanical hands can do effective massage to human cervical vertebra thoracic vertebra and lumbar vertebra.

USB-CHARGER

USB While you relax, the integrated USB port can supply, for example, your phone with power.



With one touch the quick start buttons in your armrest allow you to control the most important functions.


ARTHROSE TRACTION

Arthrosis Traction consists of gentle stretching pressure. It has a stimulating effect on the cartilages, which enhances the transport of bodily fluids, nutrients absorbation and can alleviate symptoms of arthrosis.



I-SHAPF You can expect a profound full body massage from the neck down to the buttocks.

MEMORY



ADJUSTABLE FOOTREST

The footrest can be extended, this way it is adaptable to different body sizes.

ADJUSTABLE PILLOW

If necessary, the double cushion can be placed to reduce the neck and shoulder mas-sage intensity

BACK STRETCH ROLLING Targeted stretch massages activate the body



ACUPRESSURE POINTS Key acupressure points are massaged.



and you will breathe purified oxygen during your massage.



back to ensure optimal pressure throughout the massage. The scan technology detects the neight of your shoulders, that can be modified after the scan, if required.

The massage chair is equipped with Bluetooth technology and an amphitheatre sound system that reduces ambient noise and boosts your sound experience. In this way you can enjoy your massage even with your favourite meditation application or your favourite music. We recommend starting the Bluetooth speakers with the lowest volume before playing music and then gradually increasing them to your desired level.

SLIM LCD TOUCH SCREEN REMOTE CONTROL

Jsing the touch screen, you can select your desired massage comfortably and adjust all other settings. You can navigate through the menu by clicking the buttons to select Neck/ Shoulder, Lower Back, Auto, Manual, Advance or Air Massage. Optionally, you can set the massage time to a period of 5 minutes, starting from 5 minutes up to 30 minutes.

braintronics brainwave stimulation system

The human brain is a complex organ consisting of millions of little nerve cells, called neurons. Each cell constantly produces tiny electrical impulses. Every thought, every sound and every impression, everything is transmitted and communicated by those small electrical impulses. Electricity is, so to say, the language of the brain. The sum of those impulses can be detected and diagrammed, e.g. by using an EEG. Charted on a graph this electrical activity forms a wave-like pattern, kno as brainwaves.

With braintronics[®], combining audio stimulation and guided meditations with a soothing massage, these brainwaves can be stimulated in order to reach different states of relaxation in the brain. Whether you want to sustainably enhance your concentration, or just relax and pamper yourself - everything is possible with a massage chair from the Casada family.

REFLEXOLOGY Reflexology regenerates foot muscles and has a harmonizing effect on almost all organs inside your body.


O2 IONIZER Polluted air around you will be neutralized

DUAL SENSOR BODY SCAN

The massage chair performs a body scan to capture the main areas of the neck and

















TI



CASAda

