

## SPECIAL FEATURES

**SHIATSU**  
Shiatsu (finger press) is a form of body therapy originating in Japan. The massage technique consists of soft, rhythmic, far reaching stretches and rotations. The aim of Shiatsu is to simulate individual parts of the body and to mobilise the muscles.

**3D MASSAGE**  
In this new technique, the massage robot not only moves upwards and downwards, but also forwards and backwards.

**KNEADING**  
Skin and muscles are either clasped and kneaded between the thumb and the index finger or with both hands. This massage technique is used in particular to ease tension.

**TAPPING**  
Short, tapping movements are performed with the edge of the hand, the palm or the fist. This promotes blood circulation in the skin and softens tense muscles. If tapping is performed at lung height, this can improve the release of mucus in the lung.

**KNEAD & TAP (DUAL ACTION)**  
The techniques of kneading and tapping are combined in order to relieve tension and stimulate blood circulation at the same time.

**ROLLING**  
Pleasant alternation between relaxation and releasing the muscles – also perfectly suited for a final massage, for relaxation and recovery of the muscles.

**CARBON HEATING**  
Deeply penetrating Carbon infrared heat, blood vessels expand, the blood flow increases and the muscle tone is reduced. In addition, it carries out a highly beneficial effect on nerve cells, which are responsible for reporting pain to the brain.

**TAPPING**  
Short, tapping movements are performed with the edge of the hand, the palm or the fist. This promotes blood circulation in the skin and softens tense muscles. If tapping is performed at lung height, this can improve the release of mucus in the lung.

**KNEAD & FLAP**  
A massage combination of kneading and flapping.

**AIR COMPRESSION**  
Rising and receding airbags have a pumping motion which produces natural muscle reflexes. This is particularly beneficial in the calf and foot area.

**ZERO SPACE**  
For the lying-position the massage chair moves forward and therefore requires very little distance to the wall behind. Space saving and comfortable!

**ZERO GRAVITY**  
You put yourself in the position of an astronaut – that “weightless feeling”.

**AUTOMATIC**  
Specifically developed automatic programmes for the whole body, which offer a big range of different massage techniques. Additionally you can also turn on the heating function, if you'd like.

**BACK AREA**  
An individual applicable to any size massage area.

**ADJUSTABLE FOOTREST**  
The footrest can be extended, this way it is adaptable to different body sizes.

**RECLINING POSITION**  
The adjustability of the backrest and footrest contributes to personal well-being of the body position.

**ADJUST SHOULDER**  
The massage heads can be moved upward or downward to a suitable position in the shoulder area.

**AIR INTENSITY**  
Adjustable airbag massage intensity.

**AIR AREA**  
Whole body airbag massage. Different areas are individually selectable.

**SPOT**  
The massage heads can massage a particular spot.

**SPEED**  
Massage speed is adjustable.

**WIDTH**  
The distance between the bath massage heads can be set (narrow, medium and wide).

**CALF KNEADING**  
Kneading and pressure massage functions in the leg airbags ensure that any tension in the leg muscles that may be caused by long periods of standing are relieved.

**BLUETOOTH**  
A wireless connection to your smartphone or tablet ensures a facilitated music experience.

**HUMAN HANDS FEELING**  
Specially designed massage heads imitate the hand motion of a professional masseur. The maximum stretch space is 6.5 cm, and the maximum stretch angle is 41 degrees, so the mechanical hands can do effective massage to human cervical vertebrae thoracic vertebrae and lumbar vertebra.

**SHOULDER GRASP**  
Targeted kneading massage in the shoulder/neck area – ideal for easing tension and pain in this area.

**WAIST STRETCH**  
Airbags in the lower area of the back ensure gentle stretching of the pelvis. Effectively counteracts tensions.

**VOLUME**  
Adjusting the volume.

**MEMORY**  
Using the intelligent memory function the massage chair remembers your favourite massage exactly and begins immediately without a scan.

**MUSIC**  
Further relaxation and fun with the music function.

**ADJUSTABLE PILLOW**  
If necessary, the double cushion can be placed to reduce the neck and shoulder massage intensity.

**BACK STRETCH ROLLING**  
Targeted stretch massages activate the body and have an invigorating effect on the whole body.

**ACUPRESSURE POINTS**  
Key acupressure points are massaged.

**REFLEXOLOGY**  
Reflexology regenerates foot muscles and has a harmonizing effect on almost all organs inside your body.

**HEATED ROLLERS**  
The massage chair comes with the most innovative comfortable back heating on the market so far. The rollers heat up by themselves and simulate hot stones to warm the muscles. Compared to conventional heating systems, this provides a greater surface area to apply heat from your neck all the way to the back of your thighs.

**DUAL SENSOR BODY SCAN**  
The massage chair performs a body scan to capture the main areas of the neck and back to ensure optimal pressure throughout the massage. The scan technology detects the height of your shoulders, that can be modified after the scan, if required.

**BLUETOOTH AMPHITHEATRE SOUND SPEAKERS**  
The massage chair is equipped with Bluetooth technology and an amphitheatre sound system that reduces ambient noise and boosts your sound experience. In this way you can enjoy your massage even with your favourite meditation application or your favourite music. We recommend starting the Bluetooth speakers with the lowest volume before playing music and then gradually increasing them to your desired level.

**XD FUNCTION**  
The XD massage function allows for individual adjustment of both the intensity and the spacing of the massage rollers. This ensures a particularly powerful and targeted massage – even under higher pressure or strain. Ideal for anyone seeking a deep and customizable massage experience.

**GRAPHENE WARMTH IN THE LEG AREA**  
This switches on with the heating function in the lumbar area. The Graphene technology heats the leg area within a few minutes.



**MOVABLE SIDE PARTS**  
The side parts can move forwards and backwards to provide an optimum stretching experience.

**INDIVIDUAL LEG MASSAGE**  
Separate units adapt to each leg and enable flexible treatment.

**ADAPTABLE HEAT**  
Customised heat adjustment in various areas

**175° BACKREST POSITION**  
The backrest is electrically and continuously adjustable with a tilt angle of up to 175°.

**FLAP**  
Circulation in the muscles is stimulated and the tendons are encouraged to stretch by means of light tapping with a cupped hand, which also reduces tension at the same time.



Technical Specifications	
Product Name:	Dreamliner
Voltage:	220-240 V~, 50/60 Hz
Rated Power:	240 W
Running Time:	10-20 - 30 minutes
Airbags:	Shoulders · Arms · Buttocks · Thighs · Top of Feet · Heels
Rollers:	Foot rollers
Massage Mechanism Range:	Neck to Buttocks
Remote Control:	480 X 800
Speakers:	Full-range speakers · 40 W
Weight:	167 kg
Upright Dimensions:	158 X 80.8 X 125 cm
Reclined Dimensions:	202 X 80.8 X 80.8 cm
Certificates:	CE EN13432

**Casada International GmbH**  
Obermeiers Feld 3  
33104 Paderborn

Phone: 05254 / 930980  
E-Mail: info@casada.com  
Web: www.casada.com



# DREAMLINER



# SPECIAL FEATURES

**SHIATSU**  
Shiatsu (finger press) is a form of body therapy originating in Japan. The massage technique consists of soft, rhythmic, far reaching stretches and rotations. The aim of Shiatsu is to simulate individual parts of the body and to mobilise the muscles.

**3D MASSAGE**  
In this new technique, the massage robot not only moves upwards and downwards, but also forwards and backwards.

**KNEADING**  
Skin and muscles are either clasped and kneaded between the thumb and the index finger or with both hands. This massage technique is used in particular to ease tension.

**TAPPING**  
Short, tapping movements are performed with the edge of the hand, the palm or the fist. This promotes blood circulation in the skin and softens tense muscles. If tapping is performed at lung height, this can improve the release of mucus in the lung.

**KNEAD & TAP (DUAL ACTION)**  
The techniques of kneading and tapping are combined in order to relieve tension and stimulate blood circulation at the same time.

**ROLLING**  
Pleasant alternation between relaxation and releasing the muscles – also perfectly suited for a final massage, for relaxation and recovery of the muscles.

**CARBON HEATING**  
Deeply penetrating Carbon infrared heat, blood vessels expand, the blood flow increases and the muscle tone is reduced. In addition, it carries out a highly beneficial effect on nerve cells, which are responsible for reporting pain to the brain.

**TAPPING**  
Short, tapping movements are performed with the edge of the hand, the palm or the fist. This promotes blood circulation in the skin and softens tense muscles. If tapping is performed at lung height, this can improve the release of mucus in the lung.

**KNEAD & FLAP**  
A massage combination of kneading and flapping.

**AIR COMPRESSION**  
Rising and receding airbags have a pumping motion which produces natural muscle reflexes. This is particularly beneficial in the calf and foot area.

**ZERO SPACE**  
For the lying-position the massage chair moves forward and therefore requires very little distance to the wall behind. Space saving and comfortable!

**ZERO GRAVITY**  
You put yourself in the position of an astronaut – that “weightless feeling”.

**AUTOMATIC**  
Specifically developed automatic programmes for the whole body, which offer a big range of different massage techniques. Additionally you can also turn on the heating function, if you'd like.

**BACK AREA**  
An individual applicable to any size massage area.

**ADJUSTABLE FOOTREST**  
The footrest can be extended, this way it is adaptable to different body sizes.

**RECLINING POSITION**  
The adjustability of the backrest and footrest contributes to personal well-being of the body position.

**ADJUST SHOULDER**  
The massage heads can be moved upward or downward to a suitable position in the shoulder area.

**AIR INTENSITY**  
Adjustable airbag massage intensity.

**AIR AREA**  
Whole body airbag massage. Different areas are individually selectable.

**SPOT**  
The massage heads can massage a particular spot.

**SPEED**  
Massage speed is adjustable.

**WIDTH**  
The distance between the bath massage heads can be set (narrow, medium and wide).

**CALF KNEADING**  
Kneading and pressure massage functions in the leg airbags ensure that any tension in the leg muscles that may be caused by long periods of standing are relieved.

**BLUETOOTH**  
A wireless connection to your smartphone or tablet ensures a facilitated music experience.

**HUMAN HANDS FEELING**  
Specially designed massage heads imitate the hand motion of a professional masseur. The maximum stretch space is 6.5 cm, and the maximum stretch angle is 41 degrees, so the mechanical hands can do effective massage to human cervical vertebrae thoracic vertebrae and lumbar vertebra.

**SHOULDER GRASP**  
Targeted kneading massage in the shoulder/neck area – ideal for easing tension and pain in this area.

**WAIST STRETCH**  
Airbags in the lower area of the back ensure gentle stretching of the pelvis. Effectively counteracts tensions.

**VOLUME**  
Adjusting the volume.

**MEMORY**  
Using the intelligent memory function the massage chair remembers your favourite massage exactly and begins immediately without a scan.

**MUSIC**  
Further relaxation and fun with the music function.

**ADJUSTABLE PILLOW**  
If necessary, the double cushion can be placed to reduce the neck and shoulder massage intensity.

**BACK STRETCH ROLLING**  
Targeted stretch massages activate the body and have an invigorating effect on the whole body.

**ACUPRESSURE POINTS**  
Key acupressure points are massaged.

**REFLEXOLOGY**  
Reflexology regenerates foot muscles and has a harmonizing effect on almost all organs inside your body.

**HEATED ROLLERS**  
The massage chair comes with the most innovative comfortable back heating on the market so far. The rollers heat up by themselves and simulate hot stones to warm the muscles. Compared to conventional heating systems, this provides a greater surface area to apply heat from your neck all the way to the back of your thighs.

**DUAL SENSOR BODY SCAN**  
The massage chair performs a body scan to capture the main areas of the neck and back to ensure optimal pressure throughout the massage. The scan technology detects the height of your shoulders, that can be modified after the scan, if required.

**BLUETOOTH AMPHITHEATRE SOUND SPEAKERS**  
The massage chair is equipped with Bluetooth technology and an amphitheatre sound system that reduces ambient noise and boosts your sound experience. In this way you can enjoy your massage even with your favourite meditation application or your favourite music. We recommend starting the Bluetooth speakers with the lowest volume before playing music and then gradually increasing them to your desired level.

**XD FUNCTION**  
The XD massage function allows for individual adjustment of both the intensity and the spacing of the massage rollers. This ensures a particularly powerful and targeted massage – even under higher pressure or strain. Ideal for anyone seeking a deep and customizable massage experience.

**GRAPHENE WARMTH IN THE LEG AREA**  
This switches on with the heating function in the lumbar area. The Graphene technology heats the leg area within a few minutes.



**MOVABLE SIDE PARTS**  
The side parts can move forwards and backwards to provide an optimum stretching experience.

**INDIVIDUAL LEG MASSAGE**  
Separate units adapt to each leg and enable flexible treatment.

**ADAPTABLE HEAT**  
Customised heat adjustment in various areas

**175° BACKREST POSITION**  
The backrest is electrically and continuously adjustable with a tilt angle of up to 175°.

**FLAP**  
Circulation in the muscles is stimulated and the tendons are encouraged to stretch by means of light tapping with a cupped hand, which also reduces tension at the same time.



Technical Specifications	
Product Name:	Dreamliner
Voltage:	220-240 V~, 50/60 Hz
Rated Power:	240 W
Running Time:	10-20 · 30 minutes
Airbags:	Shoulders · Arms · Buttocks · Thighs · Top of Feet · Heels · Foot rollers
Rollers:	Massage Mechanism Range: Neck to Buttocks
Remote Control:	480 X 800
Speakers:	Full-range speakers · 40 W
Weight:	167 kg
Upright Dimensions:	158 X 80.8 X 125 cm
Reclined Dimensions:	202 X 80.8 X 80.8 cm
Certificates:	CE EN 13432-1, 2014/68/EU

**Casada International GmbH**  
Obermeiers Feld 3  
33104 Paderborn

Phone: 05254 / 930980  
E-Mail: info@casada.com  
Web: www.casada.com



# DREAMLINER



A  
N  
N  
I  
V  
E  
R  
Y  
C  
A  
R  
S  
A  
R  
Y